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Ziegelei 1 D-72336 Balingen E-Mail: info@kernsohn.com Tel: +49-[0]7433- 9933-0 Fax: +49-[0]7433-9933-149 Internet: www.kernsohn.com

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KERN KXC-TM

Version 1.1 2024-11

Operating Instructions Display units

Content

1		Тес	chnical data	.4
2		Dec	claration of conformity	. 5
3		Арр	pliance overview	.6
	3.	1	Components	.6
	3.	2	Operating elements	.7
		3.2.	.1 Keyboard overview	. 7
		3.2.	.2 Numerical input	. 8
		3.2.	.3 Overview of displays	. 9
4		Bas	sic Information (General)	10
	4.	1	Proper use	10
	4.	2	Improper Use	10
	4.	3	Warranty	10
	4.	4	Monitoring of Test Resources	11
5		Bas	sic Safety Precautions	11
	5.	1	Pay attention to the instructions in the Operation Manual	11
	5.	2	Personnel training	11
	5.	3	Electrostatic sensitive components	11
6		Trai	insport and storage	12
	6.	1	Testing upon acceptance	12
	6.	1	Packaging / return transport	12
7		Unp	packing, Installation and Commissioning	13
	7.	1	Installation Site, Location of Use	13
	7.	2	Unpacking and checking	14
	7.	3	Mains connection	14
	7.	4	Rechargeable battery operation (optional)	15
		7.4.	.1 Load rechargeable battery	15
	7.	5	Connection of peripheral devices	16
	7.	6	Initial Commissioning	16
	7.	7	Adjustment of non-calibratable devices	16

7.7.1		External adjustment < ⊏用LEHE >	17
	7.7.2	External adjustment with user-defined adjustment weight < $\Box RLEud$ >	18
	7.7.3	Gravitational constant adjustment location < 다 머머리니 >	20
	7.7.4	Gravitational constant place of location < ն⊢ Rս հE >	21
7	.8 0	Calibration of calibratable devices	22
8	Verifi	cation	25
9	Basic	Operation	27
9	.1 T	urn on/off	27
9	.2 5	Simple weighing	27
9	.3 Z	Zeroing	28
9	.4 T	aring	28
9	.5 0	Change-over button and F button (standard settings)	29
	9.5.1	Switch-over weighing unit	30
	9.5.2	Display gross weight value	31
	9.5.3	Open PRE-Tare settings	31
	9.5.4	Carry out Data-Hold function	31
10	Ор	erating concept	32
11	Ар	olication <weighing></weighing>	34
1	1.1	Application-specific settings	34
1	1.2	PRE-Tare	36
	11.2.	Take over the placed weight as PRE-TARE value	36
	11.2.2	2 Enter the known tare weight numerically	37
1	1.3	Data-Hold function	37
1	1.4	Weighing Units	38
	11.4.	Setting weighing unit	38
	11.4.2	2 Weighing with multiplication factor via the application unit <ffa></ffa>	39
	11.4.3	3 Percent weighing by application unit <%>	39
12	Ар	olication <counting></counting>	40
1	2.1	Application-specific settings	40
1	2.2	Using the application	41
	12.2.	Piece counting	41
	12.2.2	2 Target counting	44
13	Ар	olication < Checkweighing >	47
1	3.1	Application-specific settings	47
1	3.2	Using the application	48
	13.2.	I Target weighing	48
	13.2.2	2 Checkweighing	51

14	Mei	nu	53
14	.1	Navigation in the menu	53
14	.2	Application menu	53
14	.3	Setup menu	54
	14.3.1	Overview < ۲۲	54
15	Cor	nmunication with peripheral devices	65
15	5.1	KERN Communications Protocol (KERN Interface Protocol)	65
15	.2	KERN alibi memory	66
15	.3	Issue functions	67
	15.3.1	Add-up mode < ╘uП >	67
	15.3.2	Data output after pressing the PRINT button < $\Pi \Pi \square \Pi \square \square \square$	69
	15.3.3	Automatic data output < 日山とロ>	70
	15.3.4	Continuous data output <	70
15	.4	Data format	71
16	Ser	vicing, maintenance, disposal	72
16	5.1	Cleaning	72
16	5.2	Servicing, maintenance	72
16	5.3	Disposal	72
17	Inst	ant help for troubleshooting	73
18	Erro	or messages	74

1 Technical data

KERN	КХС-ТМ	
Item no./ Type	TKXC-TM-A	
Display	LCD 6 digits, height 48 mm with back lighting	
Resolution (verifiable)	Single (Max.) 3000 e	
	Multi Range/Multi Intervall (Max.) 2x3000 e	
Resolution (non-verifiable)	100 – 999.999 d	
Verification class	III	
Weighing ranges	2	
Divisions	1,2,5,10, n	
DMS weighing cells	87-1100 Ω. (minimum/maximum resistance)	
Applications	Weighing, Counting, Checkweighing	
Weighing Units	g, kg, lb, pcs, %, FFA	
Allowable ambient tempera- ture	-10 °C + 40 °C	
Operating temperature range with the rechargeable battery	0 °C + 40 °C	
Humidity of air	max. 80% rel. (non-condensing)	
Electric Supply	Input voltage 100 ~ 240 V; 50 / 60 Hz; 0.4 A Overvoltage Category II Mains supply voltage fluctuations ±10 %	
	Optional storage battery TYKR-01-A (RC193650); 3.7 V; 3700 mAh	
Storage battery operation (optional)	Operating time 48 hrs (backlight off) Operating time 20 hrs (backlight on) Loading time approx. 8 hrs.	
Dimensions display unit	232 x 80 x 150 (W x D x H) [mm]	
Net weight (kg)	2,5	
Interfaces	RS-232, USB-Device, WiFi, Analogue (0-10V, 4-20mA), Ether- net, Bluetooth via KUP (optional)	
Altitude	Below 2000 m	
IP protection	IP 68	
Pollution degree	2	

* Smallest part weight for piece counting - under lab conditions:

- > There are ideal ambient conditions for high-resolution counting
- > The parts to be counted are not scattered

** Smallest part weight for piece counting - under normal conditions:

- > There are unsteady ambient conditions (draft, vibrations)
- > The parts to be counted are being scattered

2 Declaration of conformity

The current EC/EU Conformity declaration can be found online in:

www.kern-sohn.com/ce

3 Appliance overview

3.1 Components

Pos.	Designation		
1	Display		
2	Keyboard		

3.2 Operating elements

3.2.1 Keyboard overview

Button	Name	Function in Operating mode	Function in Menu
ON OFF	ON/OFF-button	 Switch on/off (press button long time) Switch on/off the display back- ground illumination (press button short time) 	
F1	F1-key	 Functions key, see chap. 9.5 	
F2	F2-key	 Functions key, see chap. 9.5 	 Navigation key Menu level back Exit menu / back to weighing mode.
ڑک پ	G-key	 Change-over button, see chap. 9.5 	 > Navigation key → > Activate menu item > Confirm selection
PRINT	PRINT-button	 Calculate weighing data via interface Display increased resolution (long keystroke, only for verified scales) 	 > Navigation key ↓ > Select menu item
→0←	ZERO-key	> Zeroing	 Navigation key ↑ Select menu item
TARE	TARE-button	➤ Taring	 Invoke application menu (press button long time)

3.2.2 Numerical input

Button	Designation	Function
		Select cipher
	Navigation key →	Confirm entry. Press button repeatedly for every digit. Wait until the numeric input window extin-guishes.
PRINT	Navigation key $oldsymbol{\Psi}$	Reduce flashing cipher (0 – 9)
→0← ↑	Navigation key ↑	Increase flashing cipher (0 – 9)

3.2.3 Overview of displays 15 14 13 12 11 10 98 M+ W1W2 PTARE HOLD PEAKNETG AUTOS MENU CAL AP - 7 0 FFA 1-2pcs% >0< 6 3lb ozkg 4-0 - 5

Position	Display	Description
1		Stability display
2	>0<	Zero display
3		Minus display
4		Tolerance marks for check weighing
5		Rechargeable battery charge indicator
6	Units display / Pcs/ %	options g, kg, lb, gn, oz or Application icon [Pcs] for piece counting or [%] for determination of percentage
7	\sim	Data transfer running
8	(î-	WIFI-symbol
9	AP	Autoprint enabled
10	Σ	Weighing data can be found in the sum memory
11	G	Display gross weight value
12	NET	Display net weight value
13	HOLD	Hold/ animal weighing function
14	PTARE	Pre-Tare function
15	W1W2	Display of selected weighing range

4 Basic Information (General)

4.1 Proper use

The balance you purchased is intended to determine the weighing value of material to be weighed. It is intended to be used as a "non-automatic balance", i.e. the material to be weighed is manually and carefully placed in the centre of the weighing plate. As soon as a stable weighing value is reached, the weighing value can be read.

It can be used in indoor and outdoor. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be imparied.

4.2 Improper Use

- Our balances are non-automatic balances and not provided for use in dynamic weighing processes. However, the balances can also be used for dynamic weighing processes after verifying their individual operative range, and here especially the accuracy requirements of the application.
- Do not leave permanent load on the weighing plate. This may damage the measuring system.
- Impacts and overloading exceeding the stated maximum load (max) of the balance, minus a possibly existing tare load, must be strictly avoided. Balance may be damage by this.
- Never operate the balance in explosive environment. The serial version is not explosion protected.
- The structure of the balance may not be modified. This may lead to incorrect weighing results, safety-related faults and destruction of the balance.
- The balance may only be used according to the described conditions. Other areas of use must be released by KERN in writing.

4.3 Warranty

Warranty claims shall be voided in case:

- Our conditions in the operation manual are ignored
- The appliance is used beyond the described uses
- The appliance is modified or opened
- Mechanical damage or damage by media, liquids, natural wear and tear
- The appliance is improperly set up or incorrectly electrically connected
- The measuring system is overloaded

4.4 Monitoring of Test Resources

In the framework of quality assurance the measuring-related properties of the balance and, if applicable, the testing weight, must be checked regularly. The responsible user must define a suitable interval as well as type and scope of this test. Information is available on KERN's home page (<u>www.kern-sohn.com</u>) with regard to the monitoring of balance test substances and the test weights required for this. In KERN's accredited calibration laboratory test weights and balances may be calibrated (return to the national standard) fast and at moderate cost.

5 Basic Safety Precautions

5.1 Pay attention to the instructions in the Operation Manual

⇒ Carefully read this operation manual before setup and commissioning, even if you are already familiar with KERN balances.

5.2 Personnel training

The appliance may only be operated and maintained by trained staff.

5.3 Electrostatic sensitive components

Electrostatic discharge (ESD) can cause damage to electronic components. Damaged components do not always lead to malfunctions immediately, but sometimes only after some time.

Therefore, take precautions for ESD protection before removing hazardous components from the packaging and carrying out work in the electronics area:

- Ground yourself before touching electronic components (ESD clothing, wristband, shoes, etc.).
- Only carry out work on electronic components at suitable ESD workstations (EPA) with suitable ESD tools (antistatic mat, conductive screwdrivers, etc.).
- Only transport electronic components outside the EPA in suitable ESD packaging.
- Never remove electronic components from their packaging if they are outside the EPA.

6 Transport and storage

6.1 Testing upon acceptance

When receiving the appliance, please check packaging immediately, and the appliance itself when unpacking for possible visible damage.

6.1 Packaging / return transport

- \Rightarrow Keep all parts of the original packaging for a possibly required return.
- \Rightarrow Only use original packaging for returning.
- ⇒ Prior to dispatch disconnect all cables and remove loose/mobile parts.
- ⇒ Reattach possibly supplied transport securing devices.
- Secure all parts such as the wind screen, the weighing plate, power supply unit etc. against shifting and damage.

7 Unpacking, Installation and Commissioning

7.1 Installation Site, Location of Use

The balances are designed in a way that reliable weighing results are achieved in common conditions of use.

You will work accurately and fast, if you select the right location for your balance.

On the installation site observe the following:

- Place the balance on a firm, level surface.
- Avoid extreme heat as well as temperature fluctuation caused by installing next to a radiator or in the direct sunlight.
- Protect the balance against direct draughts due to open windows and doors.
- Avoid jarring during weighing.
- Do not expose the device to extreme dampness for longer periods of time. Non-permitted condensation (condensation of air humidity on the appliance) may occur if a cold appliance is taken to a considerably warmer environment. In this case, acclimatize the disconnected appliance for ca. 2 hours at room temperature.
- Avoid static charge of goods to be weighed or weighing container.
- Do not operate in areas with hazard of explosive material or in potentially explosive atmospheres due to materials such as gasses, steams, mists or dusts.
- Keep away chemicals (such as liquids or gasses), which could attack and damage the balance inside or from outside.
- In the event of the occurrence of electromagnetic fields, static charges (e.g., when weighing / counting plastic parts) and unstable power supply, large display deviations (incorrect weighing results, as well as damage to the scale) are possible. Change location or remove source of interference.

7.2 Unpacking and checking

Remove device and accessories from packaging, remove packaging material and install the device at the planned workplace. Check if that there has been no damage and that all items of delivery scope are present.

Scope of delivery / serial accessories:

- Display Unit
- Operating instructions

7.3 Mains connection

Do not connect the scales to the power mains unless the information on the scales (sticker) matches the local mains voltage.

Important:

- Before starting your weighing balance, check the mains cable for damage.
- > Ensure that the power unit does not come into contact with liquids.
- > Ensure access to mains plug at all times.

7.4 Rechargeable battery operation (optional)

ATTENTION	The rechargeable battery and the battery match with each other.
	 Operating temperature range with the rechargeable battery: 0 °C + 40 °C
	The rechargeable battery can only be replaced by the same or by a type recommended by the manufacturer.
	⇒ The rechargeable battery is not protected against all environmental influences. If the rechargeable battery is exposed to certain environmental influences, it may set on fire or explode. Persons may be injured or material damage may occur.
	⇒ Protect the rechargeable battery against fire and heat.
	Do not bring the rechargeable battery in contact with fluids, chemical substances or salt.
	Do not expose the rechargeable battery to high pressure or microwaves.
	Under no circumstances the rechargeable batteries and the charging unit may be modified or manipulated.
	Do not use a defective, damaged or deformed rechargeable battery.
	Do not connect or short-circuit the electrical contacts of the rechargeable battery with metallic objects.
	⇒ Liquid may squirt out from a damaged rechargeable battery. If the liquid gets into contact with the skin or the eyes, the skin and the eyes may be irritated.
	⇒ Ensure the correct polarity when inserting or changing the recharge- able battery.
	⇒ The rechargeable battery operation is overridden when the mains adapter is connected.
	⇒ If the rechargeable battery starts to smell, being hot, changing the colour or being deformed, it must be immediately un- plugged from mains supply and from the balance if possible.

7.4.1 Load rechargeable battery

Before the first use, the rechargeable battery package should be charged for at least 15 hours.

To save the rechargeable battery, in menu (see chap. 14.3.1) the automatic switchoff function $< \exists u \models \Box F F >$ can be activated.

If the capacity of the rechargeable batteries is exhausted, $< L \Box \Box \Box \Box = 2$ appears in the display. Connect the power cable as soon as possible to load the rechargeable battery. Charging time until complete recharging is approx. 8 hrs.

7.5 Connection of peripheral devices

Before connecting or disconnecting of additional devices (printer, PC) to the data interface, always disconnect the balance from the power supply.

With your balance, only use accessories and peripheral devices by KERN, as they are ideally tuned to your balance.

7.6 Initial Commissioning

In order to obtain exact results with the electronic balances, your balance must have reached the operating temperature (see warming up time chap. 1). During this warming up time the balance must be connected to the power supply (mains, rechargeable accumulator or battery).

The accuracy of the balance depends on the local acceleration of gravity.

Strictly observe hints in chapter Adjustment.

7.7 Adjustment of non-calibratable devices

As the acceleration value due to gravity is not the same at every location on earth, each balance must be coordinated - in compliance with the underlying physical weighing principle - to the existing acceleration due to gravity at its place of location (only if the balance has not already been adjusted to the location in the factory). This adjustment process must be carried out for the first commissioning, after each change of location as well as in case of fluctuating environment temperature. To receive accurate measuring values it is also recommended to adjust the balance periodically in weighing operation.

- Prepare the required adjustment weight, see chap. 1. The adjustment weight to be applied depends on the capacity of a weighing scale. Carry out adjustment as closely as possible to admissible maximum load of weighing scale. Info about test weights can be found on the Internet at: http://www.kern-sohn.com.
 - Observe stable environmental conditions. A warm up time (see chapter 1) is required for stabilization.
 - Ensure that there are no objects on the weighing plate.
 - Avoid vibration and air flow.
 - Always carry out adjustment with the standard weighing plate in place.

The adjustment is locked in weighing scales with type approval certificate.

To disable the access lock, destroy the seal mark and actuate the adjustment switch. Position of the adjustment switch, see chap. 0. Adjustment of calibratable devices see Chap. 7.8

• Attention:

After destruction of the seal the balance must be re-verified by an authorised agency and a new verification wire/seal mark fitted before it can be reused for applications subject to verification.

7.7.1 External adjustment < c ALEHE >

⇒ Place the adjustment weight and confirm with → button, $< \exists \exists \exists \exists t \geq t$ followed by $< \neg \exists \exists \exists d \geq t$ will be displayed.

⇒ Once < ⊢ E ∩ UL d> is displayed, remove the adjustment weight.

⇒ After successful adjustment the balance automatically returns to weighing mode.
 In case of an adjustment error (e.g. objects on the weighing plate) the display will show the error message < └┌ □ ∩ └
 >. Switch off balance and repeat the adjustment process.

7.7.2 External adjustment with user-defined adjustment weight < c ALE ud >

 \Rightarrow Acknowledge selection by \rightarrow button. < $\Box E \sqcap \Box$ >,

< Put Ld > followed by the weight value of the adjustment weight to be placed will be displayed.

- ⇒ Place the adjustment weight and confirm with → button, $< \exists \exists \exists \exists b > followed by < \neg \exists \exists \exists d b > will be displayed.$
- ⇒ Once < ⊢ E ∩ UL d> is displayed, remove the adjustment weight.
- ⇒ After successful adjustment the balance automatically returns to weighing mode.

In case of an adjustment error (e.g. objects on the weighing plate) the display will show the error message $< \exists \neg \Box \neg \Box >$. Switch off balance and repeat the adjustment process.

7.7.3 Gravitational constant adjustment location < G ⊂ AAd J >

INFORMATION

- Only enter the gravitational constants after adjustment and linearisation. The two constants must be known for this.
- The two gravitational constants < **Cr ARdu**> and < **Cr Au E**> are reset to the default value after readjustment.

Set the gravitational constant at the adjustment point:

The calibration location is the location where the scale is calibrated and linearised during configuration. Before setting, find out which value of the constant is valid for you at the place of adjustment and linearisation.

INFORMATION

- Only enter the gravitational constants after adjustment and linearisation. The two constants must be known for this.
- The two gravitational constants < **[- AAdu**> and < **[- AubE**> are reset to the default value after readjustment.

Set the gravitational constant at the installation site:

The installation location is the place where the scales will be used. This enables accurate measurements. Find out which value of the constant is valid for the user before setting the scale.

7.8 Calibration of calibratable devices

INFORMATION

I

Please note that to configure a calibrated device, the calibration seal must be destroyed and the scales must be recalibrated and sealed by an authorised body (e.g. when converting to another platform).

Electric shock due to contact with live components

Electric shock leads to serious injury or death

 \Rightarrow Do not touch any live components, only the adjustment switch

NOTE

⇒ Please be sure to observe the notes on electrostatically sensitive components in the chapter "Electrostatic sensitive components".

Open the display unit:

1. Loosen the screws on the back of the display unit.

2.

⇒ Make sure that you do not damage any cables (e.g. by tearing them off or pinching them).

NOTE

Carefully open both halves of the display unit.

Open service menu:

ON OFF

+

circuit board (for the position of the adjustment switch, see Chap.0)
 ⇒ Switch on the appliance and press the adjustment switch

⇒ Remove the cover of the adjustment switch on the

- \Rightarrow Wait until **<H** 10 > appears on the display
- ⇒ Release buttons
- ⇒ The device can now be configured in calibrated mode

Carry out adjustment

H 1()

To carry out the adjustment or set the gravitational constants, the service menu must be opened and called up via $\langle \exists d d u \exists b \rangle \rightarrow \langle \Box \exists b \rangle$. Operation is then the same as in chapter 7.7.1 described.

Note

- ⇒ Make sure that you do not damage any cables (e.g. by tearing or pinching them).
- \Rightarrow Make sure that any seals present are in their intended place.
- **1.** Carefully fold together both halves of the display unit.
- **2.** Screw the display unit together (tightening torque = $5 \text{ Nm} \pm 5\%$).

8 Verification

General:

According to EU directive 2014/31/EU balances must be officially verified if they are used as follows (legally controlled area):

- For commercial transactions if the price of goods is determined by weighing.
- For the production of medicines in pharmacies as well as for analyses in the medical and pharmaceutical laboratory.
- For official purposes
- For manufacturing final packages

In cases of doubt, please contact your local trade in standard.

Balances in the legally controlled area (-> verified balances) must keep the error limits in the verification validity period – normally they are the double of the verification error limits.

When this verification validity period expires, a re-verification must be carried out. Should be necessary an adjustment of the balance to keep the verification error limits to satisfy the reverification requirements, this is not deemed a warranty case.

Verification notes:

An EU type approval exists for balances described in their technical data as verifiable. If the balance is used where obligation to verify exists as described above, it must be verified and re-verified at regular intervals.

Re-verification of a balance is carried out according to the respective national regulations. The validity for verification of balances in Germany is e.g. 2 years.

The legal regulation of the country where the balance is used must be observed!

Yerification of the balance is invalid without the seal.

The seal marks attached on balances with type approval point out that the balance may only be opened and serviced by trained and authorised specialist staff. If the seal mark is destroyed, verification looses its validity. Please observe all national laws and legal regulations. In Germany a re-verification will be necessary.

Position of the official seals:

9 Basic Operation

9.1 Turn on/off

Start-up:

est.

Press the **ON/OFF** button.

Wait until the weight display appears

The scales are now ready for operation using the last active application

The display lights up and the balance carries out a selft-

Switching off:

Keep **ON/OFF** button pressed until the display disappears

9.2 Simple weighing

- Check zero display [>0<] and set to zero with the help of the **ZERO** key, as required.
- Place goods to be weighed on balance
- Wait until the stability display appears (M).
- Read weighing result.

Overload warning

Overloading exceeding the stated maximum load (max) of the device, minus a possibly existing tare load, must be strictly avoided.

This could damage the instrument.

Exceeding the maximum load is indicated by the display " $\int - \frac{1}{2}$ ". Unload balance or reduce preload.

9.3 Zeroing

In order to obtain optimal weighing results, reset to zero the balance before weighing. Zeroing is only possible in the range $\pm 2\%$ Max.

For values greater than ± 2% maximum the error message < 2L $_{1}\Pi$ $_{1}E>$ is displayed

- Unload the balance
- Press the **ZERO** key to set the balance to zero.

9.4 Taring

The dead weight of any weighing container may be tared away by pressing a button, so that the following weighing procedures show the net weight of the goods to be weighed.

Put weighing container on the weighing plate.

- Wait until the stability display appears ► →), then press TARE key. The weight of the container is now internally saved. Zero display and indicator <NET> will appear. <NET> informs that all shown weight values are net values.
- When the balance is unloaded the saved taring value is displayed with negative sign.
 - To delete the stored tare value, unload the weighing plate and press the **TARE** key or the **ZERO** key.
 - The taring process can be repeated any number of times, e.g. when adding several components for a mixture (adding). The limit is reached when the taring range capacity is full.
 - Numerical input of tare (PRE-TARE)

9.5 Change-over button and F button (standard settings)

The change-over button Rand the F button can be allocated with different functions.

The following functions are set as per standard (< dEFAuL>) in the different weighing applications:

R ₽	Short key pressing	Long key pressing
HE IH	 When pressed for first time: Setting weighing unit Switch-over between the weighing units 	Display gross weight value
count	 When pressed for first time: Setting the reference quantity Switch-over between the weighing units 	When the balance has been tared and the weighing unit is displayed, you can change the display be- tween gross weight, net weight and tare weight by pressing the button long time.
chEcĥ	 When pressed for first time: Setting weighing unit Switch-over between the weighing units 	When the balance has been tared and the weighing unit is displayed, you can change the display be- tween gross weight, net weight and tare weight by pressing the button long time.

F1	Short key pressing	Long key pressing
HE ih	Open PRE-TARE settings	 Carry out Data-Hold function
count	 Setting the reference quantity 	No function assigned
chEcĥ	 Open settings for check- weighing 	 Open settings for target weighing

F2	Short key pressing	Long key pressing
비도 나		
count	 Select application 	No function assigned
chEch		

The standard settings ($< dEFR_{u}L >$) for the <Weighing> application are described below.

9.5.1 Switch-over weighing unit

As per standard the change-over button \widehat{R} is set so that is it possible to switch-over between the weighing units by **shortly** pressing.

Enable unit:

The unit for quick selection can be determined when the $rac{2}$ -button is shortly pressed for the first time.

- ⇒ Press the button and wait until the display flashes.
- ⇒ Use the navigation keys ↓↑ to select the weighing unit and confirm on → button.

Switch over unit:

Using $\overrightarrow{\sim}$ button, it is possible to switch over between the enabled unit 1 and unit 2.
Enable another unit:

- Select menu setting < u□ i ≥ and confirm on → button.
- \Rightarrow Wait until the display flashes.
- ⇒ Use the navigation keys ↓↑ to select the weighing unit and confirm on → button.

For the required settings of an application unit (%, FFA) selection, please see chap. 11.4.2 and 11.4.3.

9.5.2 Display gross weight value

As per standard the change-over button $\widehat{\sim}$ is set so that is it possible to display the gross weight value by long-time pressing.

⇒ Keep the the button pressed until the display shows the gross weight value.

After releasing the button, the gross weight value will be kept in the display for a short time.

9.5.3 Open PRE-Tare settings

As per standard the F-key is set so that the menu setting $< P \vdash \exists \neg E >$ is invoked by pressing the key **shortly**. Further settings, see chap. 11.2.

9.5.4 Carry out Data-Hold function

As per standard the F-key is set so that the Data-Hold function $< h \Box L d >$ is executed by pressing the key for a long time, see chap. 11.3.

10 Operating concept

From factory the balance is delivered with various applications (weighing, check weighing, counting). After the first start-up the balance is in the <Weighing> application.

In the **application menu** (see chap. 14.2) however, you can define, selecting an application, in which mode the balance after switching-on has to continue working. Either as per standard in weighing mode or e.g. in check mode or counting mode.

Selecting an application:

According to the selected application in the application menu just appear the application-specific settings, so that you reach the target quickly without deviation.

- Information about the application-specific settings you will find in the description of the respective application.
 - All basic settings and parameters, which influence the whole operation of the balance, are resumed in the **Setup Menu** (see chap. 14.3.1) These settings remain valid for all applications.
 - The number of the available applications depends on the model.

Change application:

- ⇒ Press the TARE button and keep it pressed until the first menu item of the setup menu will be displayed
- ⇒ Use the \checkmark button to select the menu setting < $\Pi \Box dE$ > and acknowledge with → button. The current setting will be displayed.
- \Rightarrow Press the \checkmark button to select the required unit and confirm by pressing the \rightarrow button.

11 Application < Weighing>

How to carry out a simple weighing and taring, please refer to chap. 9.2 or 9.4. Further specific settings you will find in the following chapters.

Shouldn't the application <Weighing> already be enabled, select the menu setting < $\Pi \Box dE > \Rightarrow < \exists E \ (h >)$, see chap. 10

11.1 Application-specific settings

Call up menu:

- \Rightarrow Press the **TARE** key and hold it until < $P \subseteq \Pi \in \square >$ is displayed.
- \Rightarrow The display changes to < $\square \square \square \square$ followed by < $\square \square \square$ = .
- ⇒ Navigation in menu see chap. 14.1

Overview (not verifiable models):

Level 1	Level 2	Level 3	Description / Chapter		
	ActuAL	Take over the placed weight as PRE-TARE value,, see chap. 11.2.1			
PRE-TARE	ՈጸոսጸԼ	Numerical input of the tare weight, see chap. 11.2.2			
	cLEAr	Delete PRE-TARE value			
hold	-	Start-Hold fur	Start-Hold function, see chap. 11.3		
שה ו ב Units	available weigh- ing units, see chap. 1	igh- eeThis function defines in which weighing unit the result will b displayed, see chap. 11.4.1Application unit counting			
	pcs				
	FFA	Multiplication factor see chap. 11.4.2			
	%	Application unit for determining percentages see chap. 11.4.3			
NodE	HE 'H	Weighing			
Applications	count	Counting see chap. 10			
	chEch	Check weighing			

Overview (verifiable models):

Level 1	Level 2	Level 3	Description / Chapter	,
	ActuAL	Take over the placed weight as PRE-TARE value, see 11.2.1Numerical input of the tare weight, see chap. 11.2.2		E value, see chap.
PRE-TARE	NANUAL			ар. 11.2.2
	cLEAr	Delete PRE-TARE value		
hold	-	Start-Hold function, see chap. 11.3		
טה וב	g	This function defines in which weighing unit the result will I		
units	kg	displayed, see chap. 11.4.1		
NodE	HE 'P	Weighing		
Applications	count	Counting see chap. 10		
	chEch	Check weighing		

11.2.1 Take over the placed weight as PRE-TARE value

< PEArE > = < ActuAL >

The entered tare weight remains valid until a new tare weight is input. To delete press the TARE key or confirm the menu setting $< \Box L \Box \Box = 0$ using the \rightarrow button.

11.2.2 Enter the known tare weight numerically

 $< PEArE > \rightarrow < NAnuAL >$

The entered tare weight remains valid until a new tare weight is input. To delete enter the zero value or confirm the menu setting $< \Box L \Box \Box = 0$ using the \rightarrow button.

11.3 Data-Hold function

- ⇒ Menu setting < hoLd >
- \Rightarrow Place goods to be weighed.
- \Rightarrow Acknowledge by \rightarrow button.

⇒ The first stable weight value is symbolised by [HOLD] in the upper edge of the display. After the load is removed, the value is left in the display for another 10 seconds.

11.4 Weighing Units

11.4.1 Setting weighing unit

1

- Select menu setting < u□ ↓ ≥ and confirm on → button.
- \Rightarrow Wait until the display flashes.
- ⇒ Use the navigation keys ↓1 to select the weighing unit and confirm on → button.

- For the required settings of an application unit (FFA, %) selection, please see chap. 11.4.2 and 11.4.3.
 - Using the
 button (standard setting) you can switch between the active unit 1 and unit 2 (standard setting of buttons, see chap. 9.5. Other setting options, see chap. 14.3.1.

11.4.2 Weighing with multiplication factor via the application unit <FFA>

Here you determine with which factor the weighing result (in gram) will be multiplied.

By that way, e.g. a known error factor in the weight determination can be immediately taken into account.

11.4.3 Percent weighing by application unit <%>

The application unit <%> allows to check the weight of a sample in percent, based on a reference weight.

12 Application <Counting>

1

Shouldn't the application <Counting> already be enabled, select the menu setting < $\Pi \Box dE > \Rightarrow < \Box \Box \Box \Box E >$, see chap. 10

12.1 Application-specific settings

Call up menu:

- \Rightarrow Press the **TARE** key and hold it until < $PP \subseteq \Pi E \Box$ > is displayed.
- \Rightarrow The display changes to < $\Box \Box \Box \Box \Box \Box d$ > followed by < $\Box EF$ >.
- ⇒ Navigation in menu see chap. 14.1

Overview:

Level 1	Level 2	Level 3	Description / Ch	napter			
r E F	5	Reference quantity 5					
Reference quantity	10	Reference quantity	Reference quantity 10				
	20	Reference quantity 20					
	50	Reference quantity 50					
	FrEE	Optional, numerical input, see chap. 3.2.2					
	տԲսե	Input of piece weight, numerical input, see chap. 3.2.2					
PER-E PRE-TARE	ActuAL	Take over the placed weight as PRE-TARE value, see chap.11.2.1Numerical input of the tare weight, see chap.11.2.2					
	∩A∽∪AL						
	cLEAr	Delete PRE-TARE value					
£ArGE£	UALUE	Counting mode					
Target counting	ΕΓΓυΡΡ	Upper tolerance		see chan 1222			
	ErrLoU	Lower tolerance					
	cLEAr	Delete settings					
NodE	count Counting						
Applications	chEch	Check weighing see chap. 10					
	HE 'P	Weighing					

12.2 Using the application

12.2.1 Piece counting

Before the balance can count parts, it must know the average part weight (i.e. reference). Proceed by putting on a certain number of the parts to be counted. The balance determines the total weight and divides it by the number of parts, the so-called reference quantity. Counting is then carried out on the basis of the calculated average piece weight.

- The higher the reference quantity the higher the counting exactness.
 - Especially high reference must be selected for small parts or parts with considerably different sizes.
 - Smallest counting weight see table "Technical data".

1. Set reference

Reference quantity 5, 10, 20 or 50:

- FEF

 □

 □

 □

 □

 □

 □

 □

 □

 □

 □

 □

 □

 □

 □

 □

 □

 □

 </tabl
 - \Rightarrow If necessary, put on and tare the weighing container.
 - \Rightarrow Put on the desired quantity of reference pieces.
 - \Rightarrow Invoke menu setting < $\neg EF$ > and confirm by \rightarrow button.
 - Use the navigation keys ↓1 to select the reference piece quantity (5, 10, 20, 50) according to the placed reference and confirm with the → button.
 - The balance will calculate the average item weight and then displays the quantity of pieces.
 - Remove reference weight. The balance is now in piece counting mode counting all units on the weighing plate.

Reference quantity user-defined:

\Rightarrow If necessary, put on and tare the weighing container. -EF \Rightarrow Put on the desired quantity of reference pieces. \Rightarrow Invoke menu setting < ref > and confirm by \rightarrow button. pcs \Rightarrow Use the navigation keys \downarrow to select the setting < $F \neg EE$ > and confirm on \rightarrow button. \Rightarrow The numeric input window appears. \Rightarrow Enter and confirm the quantity of the placed reference parts, numerical input see chap. 3.2.2 00033 ⇒ The balance will calculate the average item weight and pcs then displays the quantity of parts. $\mathbf{\Psi}$ ⇒ Remove reference weight. The balance is now in piece pcs pcs counting mode counting all units on the weighing plate.

Counting with optional piece weight:

- ⇒ Invoke menu setting < $\neg EF$ > and confirm on → button.
- Solution Solution Solution Solution Solution Solution → Use the navigation keys ↓1 to select the weighing unit and confirm on → button.
- ⇒ Use the navigation keys ↓1 to select the comma position and confirm on → button.

- ⇒ Enter piece weight, numerical input see chap. 3.2.2, the active digit flashes.
- \Rightarrow Acknowledge by \rightarrow button.

The balance is now in piece counting mode counting all units on the weighing plate.

2. Parts counting

 \Rightarrow If necessary, put on and tare the weighing container.

⇒ Fill the counting quantity. The piece quantity is shown directly in the display.

1 Use the *i* to switch between piece quantity and weight display (standard setting see chap. 9.5).

12.2.2 Target counting

The <Target counting> application variant allows weighing of goods within set tolerance limits in keeping with a determined target quantity.

Reaching the target quantity is indicated by an acoustic (if activated in menu) and an optic signal (tolerance marks).

Optical signal:

The tolerance marks provide the following information:

♠	Target quantity exceeds defined tolerance
OK	Target quantity within defined tolerance
LO.	Target quantity below defined tolerance

Acoustic signal:

The acoustic signal depends on the menu setting $< \Box E \sqcup P \Rightarrow \Box E E P E \sqsubset >$, see chap. 14.3.1.

Procedure:

1. Define target quantity and tolerances

ErrLoy

- ⇒ Use the navigation keys \downarrow 1 to select the setting < $E \neg \neg \neg$ L \Box H> and confirm on → button.
- ⇒ Use the navigation keys ↓1 to select the weighing unit and confirm on → button.
- The numeric input window appears. The active digit is flashing.
- ⇒ Enter the lower tolerance (for numerical input, see chap. 0) and confirm the entry.
- ⇒ The balance returns to the < ErrLo \square > menu.

 \Rightarrow Press repeatedly \leftarrow button to exit menu.

Finished the setting works, the weighing balance will be ready for target counting.

2. Start tolerance check:

- ⇒ Determine the average piece weight, see chap. 12.2.1
- ⇒ Place the weighed material and check by means of the tolerance marks / acoustic signal if the weighed material is within the defined tolerance.

1

The entered values will remain valid until new values are entered.

13 Application < Checkweighing >

Shouldn't the application <Checkweighing> already be enabled, select the menu setting < $\Pi \Box dE > \Rightarrow < \Box h E \Box h >$, see chap. 10

13.1 Application-specific settings

Call up menu:

- ⇒ Press the **TARE** key and hold it until < P ⊂ Π ⊆ ∩ > is displayed.
- The display changes to $< \Box h \exists \Box d >$ followed by $< L : \Box : \Box : b >$.
- ⇒ Navigation in menu see chap. 14.1

Level 1	Level 2	Level 3	Description / Ch	apter	
£ArGE£	UALUE	Target weight, numerical input, see chap. 3.2.2			
Target weighing,	ΕΓΓυΡΡ	Upper tolerance, numerical input see chap. 3.2.2			
see chap. 13.2.1	ErrLoU	Lower tolerance, numerical input see chap. 3.2.2			
	cLEAr	Delete settings			
L M LS	լ "ՈսԲԲ	Upper limit value, n	umerical input see cha	p. 3.2.2	
check weighing, see chap. 13.2.2	լ սլեօզ	Lower limit value, numerical input see chap. 3.2.2			
	cLEAr	- Delete settings			
PER-E PRE-TARE	ActuAL	Take over the placed weight as PRE-TARE value, see chap.11.2.1			
	NAnuAL	Numerical input of the tare weight, see chap. 11.2.2			
	cLEAr	Delete PRE-TARE	value		
NodE	HE 'H	Weighing			
Applications	count	Counting see chap. 10			
	chEch	Check weighing			

Overview:

13.2 Using the application

13.2.1 Target weighing

The <target weighing> application variant allows weighing of goods within set tolerance limits in keeping with a determined target weight.

Reaching the target weight is indicated by an acoustic (if activated in menu) and an optic signal (tolerance marks).

Optic signal:

The tolerance marks provide the following information:

f	Upper limit	
ок	Target weight	
LO	Lower limit	

Acoustic signal:

The acoustic signal depends on the menu setting $< \Box E \sqcup P \Rightarrow \Box E E P E \sqsubset >$, see chap. 0.

Procedure:

1. Define target weight and tolerances

3. Start tolerance check:

⇒ Place the weighed material and check by means of the tolerance marks / acoustic signal if the weighed material is within the defined tolerance.

Load below specified to- lerance	Load within specified to- lerance	Load exceeds specified tolerance	
G 0.9864 kg			

The entered values will remain valid until new values are entered.

To delete the values, select menu setting $< \Box \Box \Box \Box \Box = < \Box \Box \Box \Box = >$ and confirm on \rightarrow button.

13.2.2 Checkweighing

With the <Checkweighing> application variant you can check if the weighing good is within a predefined tolerance range.

When limit values are exceeded below or above, an acoustic signal (if enabled in menu) will sound and an optic signal (tolerance marks) will be displayed

Optic signal:

The tolerance marks provide the following information:

	Weighed-in goods exceed predefined tolerance
οκ	Weighed-in goods within predefined tolerance
LO	Weighed-in goods below predefined tolerance

Acoustic signal:

The acoustic signal depends on the menu setting $< \Box E \sqcup P > \Rightarrow < \Box E P E \sqsubset >$, see chap. 14.3.1.

Procedure:

1. Define limit values

⇒ Using the navigation keys ↓1 select the setting
 Select < L ,□ , L , ⇒ and confirm on → button.

< L INUPP > will appear.

- ⇒ Press → button to confirm, the numeric input window for entering the upper limit value will appear. The active digit is flashing.
- ⇒ Enter upper limit value (numerical input see chap. 3.2.2) and confirm the entry.

The balance returns to the < L $\Pi \cup PP >$ menu.

 \Rightarrow Use the navigation keys \downarrow to select setting < L \square L \square H >.

- ⇒ Press → button to confirm, the numeric input window for entering the lower limit value will appear. The active digit is flashing.
- ⇒ Enter lower limit value (numerical input see chap. 3.2.2) and confirm the entry.

The balance returns to the < L $\Pi L \Box H >$ menu.

 \Rightarrow Press repeatedly \leftarrow button to exit menu.

Finished the setting works, the weighing balance will be ready for checkweighing.

2. Start tolerance check:

⇒ Place the weighed material and check by means of the tolerance marks / acoustic signal if the weighed material is within the defined tolerance.

Load below specified to- lerance	Load within specified to- lerance	Load exceeds specified tolerance	
G.8864 kg			

The entered values will remain valid until new values are entered.

1

14 Menu

14.1 Navigation in the menu

Call up menu:

Select and adjust parameters:

Scrolling on one level	Use the navigation buttons to select the individual menu blocks one by one.			
	Use the navigation key $oldsymbol{\Psi}$ to scroll down.			
	Use the navigation key \clubsuit to scroll up.			
Activate menu item / Confirm selection	Press navigation key ->			
Menu level back / back to weighing mode	Press navigation key 🗲			

14.2 Application menu

1

The application menu allows you a fast and targeted access to the respectively selected application (see chap. 10).

An overview of the application-specific settings you will find in the description of the respective application.

14.3 Setup menu

In the setup menu you have the possibility to adapt the behaviour of the balance to your requirements (e.g. environmental conditions, especial weighing processes).

These settings are global and do not depend on the selected application (with exception of: <bubble un>).

14.3.1 Overview < SELuP >

Not verifiable models:

		other level	s / description	
	Level 2	Descriptio	n	
cAL	cALEHE	➔ External	adjustment, see chap. 7.7.1	
Adjustment	cALEud	→ External	adjustment, user-defined, see chap. 7.7.2	
	Cr AAdJ	→ Gravity	constant adjustment site, see chap.7.7.3	
	նոԶսեն	→ Gravity	constant installation site, see chap. 7.7.4	
соП	-5232	bAud	600	
Communication	¢		1500	
	ubb-d		2400	
			4800	
			9600	
			14400	
			19200	
			38400	
			57600	
			1 15200	
			128000	
			256000	
		4858	ባለሁ የደረጉ	
			846 .65	
		PAr ity	nonE	
			odd	
			EUEn	
		ъεοΡ	156 it	
			256 /65	
		hAndSh	nonE	
		Protoc	ĥcP	
	A_A-0P	0- 100	0 – 10 V	
		4-20NA	4 – 20 mA	
	8680	on	WLAN ein	
	oFF		WLAN aus	

	1		1		1		
Pr int	IntFeE		-5232		RS 232 interface*		
Data output	t		ubb-d		USB interface*		
			<u>ы п_</u>		WLAN interfa	ace*	
			0100	BLHU		ection with KUP interface	
	აიი		00		Switch on / o	ff add-up mode,	
			oFF		see chap. 15	.3.1	
	PrNodE	ברים			on, oFF		
			IHnuH	ΠΑσυΑΓ		by pressing the n, see chap. 15.3.2	
			RutoPr	-	on, oFF		
					Automatic da	ta output with stable and	
					positive weig	hing value	
					ter zero displ	av and stabilisation, de-	
					pending on th	ne settings	
					< 2- AnGE	>, selectable	
					(off, 1, 2, 3,4	,5). < L'r HnuE > defines	
					with d results	the factor for d. This factor multiplied	
					exceeded, a	value cannot more be con-	
					sidered as stable.		
				oFF	Continuous c	lata output	
					SPEEd	Setting output interval see chap. 15.3.4	
					2Ero	on, oFF	
			cont on			0 (unloaded) also transmit	
						continuously	
					SEAPLE	on, oFF	
						Transmit stable values only	
		8E (Ghe	56LPrt		on, oFF	Displayed weight value is transmitted	
					ნითაა	on, oFF	
					nEt	on, of F	
					EALE	on, oFF	
			GntPrt		ForNAL	LonG (detailed measu- rement protocol)	
						שלים surement protocol)	

LΑΥουΈ	nonE	on, oFF Standard lay	out
		NodEL	on, oFF
			Output model designation of the scale
		SEr iAL	on, oFF
	ubEr		Output serial number of the scale
		AL id	Alibi-ID output
		94FE	Date output
		F 'UE	Time output
		on, oFF	
		GLP complia	nt weighing protocol output
- 5 - 5 -	no	Do not delete	e settings
rtbtt	965	Delete settin	gs

BEEPE r	REYS	oFF	Switch on / off	acoustic signal by pressing
Acoustic signal		on	button	
	chEcĥ		oFF	Acoustic signal off
			5608	Slow
		ch-oh	560	Standard
			FASE	Fast
			cont.	Continuous
			oFF	Acoustic signal off
			5608	Slow
		ch-Lo	560	Standard
			FASE	Fast
			cont.	Continuous
		ch-h ,	oFF	Acoustic signal off
			5608	Slow
			560	Standard
			FASE	Fast
			cont.	Continuous
AutoFF	NodE	oFF	Automatic switch-off function switched off	
Automatic switch-off function in rechargeable battery operation		Auto	The balance is automatically switched-off according to the time without load change or without operation defined in menu item \pm $.\Pi E$ >	
		only0	Automatic swit	ch-off only with zero display
	ι ΠΕ	305	After the set tir	ne without load change or
		111 m	operation the b	palance will switch off auto-
		2 N m	matically	
		5 M in		
		60 N in		

button Key allocation			dEFAul	Standard settings, see chap. 9.5
			oFF	Button disabled
			บก เป	Set weighing unit, see chap. 11.4.1
			NodE	Select weighing application, see chap. 10
			hold	Execute HOLD function, s.Kap. 11.3 *only for the application <weighing></weighing>
	F I-REY ¢	ამიან	0-0-0	Open PRE-Tare settings, see chap. 11.2
	F2-REY ¢	¢ LPuSh	FERFE	<pre>>Weighing>, <check- weighing></check- </pre>
chAnGE		- F F	Set reference quantity, see chap. 12.2.1	
			' _'	*only for the application <counting></counting>
			լ ուեշ	Open settings for checkweighing, see chap. 13.2.2
				*only for the application <checkweighing></checkweighing>
			ЕЯгСЕЕ	Open settings for target weighing, see chap. 13.2.1 *only for the application <checkweighing></checkweighing>
bL ,БhE Display background il-	NodE	ALUAYS	Background lig on permanent	ghting of display is switched y
lumination		E MEr	The background illumination is automatically switched-off according to the time without load change or without operation defined in menu item $< E (\Pi E >$	
		no bL	Display backg switched off	lay background illumination always ched off
	F 'UE	55 105	Definition, afte illumination is	r which time the background automatically switched-off
		<u>כטכ</u> ח ח ו		ange of without operation.
		<u>50 m</u>	-	
		<u>50 m</u> 300 m	-	

ERFEFE Taring range	IDD% ¢ ID%	Definition max. taring range, selectable 10% - 100%. Numeri- cal input see chap. 3.2.2		
2trAch	on	Automatic zero tracking [<u><</u> 3d]		
Zerotracking	oFF	In the event that small quantities are removed or add to the material to be weighed, incorrect weighing re- sults can be displayed due to the "stability compensa tion". (e.g. slow flow of liquids from a container place on the balance, evaporating processes).		
			When appo it is advisat	rtioning involves small variations of weight, ole to switch off this function.
ARF 'UE	566	-20)22-	Enter the year
Date and time	dAForN	15-31		Enter the month and the day
	t For N	239	j9.59	Enter the time (hours, minutes, seconds)
սուեն Units	available weighing units / appication units, see chap. 1	ם, םFF Using this function you can define which weighing units are available in the application-specific menu < עם ול>. The units selected by < םם > are available in the application-spe- cific menu.		
NodE's	8E .h	Weig	hing	
weighing applications	count	Cour	nting	
	chEcĥ	Check weighing		
Loch	SELLCH	□n, □FF This function can be used to block access to the setup menu. When set to <□□>, the entry of a 6-digit number is required as a password.		
rESEE	Reset balance s	nce settings to factory settings		

Verifiable models:

		other leve	ls / description	
Lever		Description		
coN	-5232	bRud	600	
Communication	Ф		1200	
	ubb-d		2400	
			4800	
			9600	
			14400	
			19200	
			38400	
			57600	
			1 15200	
			158000	
			256000	
		98F8	7d6 /65	
			8d6 /65	
		PAr 169	nonE	
			odd	
			EUEn	
		SEOP	lbb ib	
			256 /65	
		hAndSh	nonE	
		Protoc	ϜϲP	
	8~8-oP	0- 100	0 – 10 V	
		4-20NA	4 – 20 mA	
	8LAv	on	WLAN ein	
		oFF	WLAN aus	

	I		1			
Pr int	IntFeE		<u>r5232</u>		RS 232 interf	ace*
Data output	Data output		სემ-950		USB interface	e*
					WLAN interfa	ace*
			BEHN		*only in conn	ection with KUP interface
	500		on		Switch on / o	ff add-up mode,
			oFF		see chap. 15	.3.1
	PrNodE	ברים			on, oFF	
			NA~uAl		Data output b PRINT buttor	by pressing the n, see chap. 15.3.2
			RutoP	-	on, oFF	
					Automatic da positive weig see chap.15. ter zero displ pending on tt < Tr RnEE (off, 1, 2, 3,4 the factor for with d results exceeded, a sidered as st	ta output with stable and hing value 3.3. Another output only af- ay and stabilisation, de- ne settings >, selectable ,5). < 고-RnGE > defines d. This factor multiplied in the threshold; when it is value cannot more be con- able.
				oFF	Continuous c	lata output
					SPEEd	Setting output interval see chap. 15.3.4
			CONC	on	CEro	ם ה, םFF 0 (unloaded) also transmit continuously
		86 _{(Сн} е	SGLPrt		on, oFF	Displayed weight value is transmitted
					նունե	on, oFF
				nEE	on, oFF	
				FUL	on.oFF	
			GntPri	E	ForNAL	LonG (detailed measu- rement protocol)
					שלם (standard mea- surement protocol)	

	LAYout	nonE	on, oFF Standard lay	out
			NodEL	on, oFF
				Output model designation of the scale
			SEr AL	on, oFF
		ubEr		Output serial number of the scale
			AL id	Alibi-ID output
			948FE	Date output
			F 'UE	Time output
			on, oFF	
			GLP complia	nt weighing protocol output
	rESEE	no	Do not delete	e settings
		YE5	Delete setting	gs

ЬЕЕРЕг	REYS	oFF	Switch on / off	acoustic signal by pressing	
Acoustic signal		on	button		
	chEch		oFF	Acoustic signal off	
			5608	Slow	
		ch-oh	5Ed	Standard	
			FASE	Fast	
			cont.	Continuous	
			oFF	Acoustic signal off	
			5608	Slow	
		ch-LO	560	Standard	
			FASE	Fast	
			cont.	Continuous	
		ch-h (oFF	Acoustic signal off	
			5608	Slow	
			560	Standard	
			FASE	Fast	
			cont.	Continuous	
AutoFF	NodE	oFF	Automatic switch-off function switched off		
Automatic switch-off function in rechargeable battery operation		Ruto	The balance is according to th or without ope L .⊓E >	automatically switched-off the time without load change ration defined in menu item <	
		onlyO	Automatic swit	Automatic switch-off only with zero display	
	F 'UE	305	After the set tir	me without load change or	
		II .m	operation the b	palance will switch off auto-	
		<u> </u>	matically		
		<u>50 m</u>	-		
		<u></u>	-		
		60 M m			

button Key allocation			dEFAul	Standard settings, see chap. 9.5	
			oFF	Button disabled	
			טה וב	Set weighing unit, see chap. 11.4.1	
			NodE	Select weighing application, see chap. 10	
			hold	Execute HOLD function, s.Kap. 11.3 *only for the application <weighing></weighing>	
	FI-REY			Open PRE-Tare settings, see chap. 11.2	
	¥ F2-KEY ¢	srush ¢ LPush	PEArE	*only for the applications <weighing>, <check- weighing></check- </weighing>	
	chAnGE		-EE	Set reference quantity, see chap. 12.2.1	
				*only for the application <counting></counting>	
			L n es	Open settings for checkweighing, see chap. 13.2.2 *only for the application	
			ЕЯгСЕЕ	Open settings for target weighing, see chap. 13.2.1 *only for the application <checkweighing></checkweighing>	
占L ・Gトヒ Display background il-	NodE	ALUAYS	Background lig on permanent	hting of display is switched y	
lumination		צ יח∃ר	The background illumination is automatically switched-off according to the time without load change or without operation defined in menu item $< E$ ($\Pi E >$		
		no bl	Display backgi switched off	splay background illumination always itched off	
		55 105	Definition, afte illumination is a	r which time the background automatically switched-off	
		305	without load ch	nange or without operation.	
	EdlE				
		<u>ciin</u> 50	-		
		300			

98F 'UE	586	-2055-	Enter the year		
Date and time	dAForN	15-31	Enter the month and the day		
	E For N	2359.59	Enter the time (hours, minutes, seconds)		
սուեն	available	on, oFF			
Units	weighing units / appication units, see chap. 1	n you can define which weighing units are oplication-specific menu < uロ ルン. The < ロロ > are available in the application-spe-			
NodES Waishing applications	8E ih	Weighing			
	count	Counting			
	chEch	Check weighing			
		on, oFF			
ιοςδ	SEELch	This function can be used to block access to the setup menu. When set to $<\Box\Box$, the entry of a 6-digit number is required as a password.			
rESEE	Reset balance s	ince settings to factory settings			

15 Communication with peripheral devices

15.1 KERN Communications Protocol (KERN Interface Protocol)

KCP is a standardized set of interface orders for KERN balances, which allows many parameters and device functions to be called up and controlled. KERN devices that have KCP can use it to connect easily to computers, industrial control systems and other digital systems. A detailed description you will find in the "KERN Communications Protocol" manual, available in the download area on our KERN homepage (www.kern-sohn.com).

To activate KCP please observe the menu overview of your balance's operating instructions.

KCP is based on simple ASCII orders and replies. Every interaction consists of an order, possibly with arguments separated by spaces and finished by <CR>< LF>.

The KCP orders supported by your balance may be queried emitting the order "I0" followed by CR LF.

10	Shows all implemented KCP orders
S	Sending stable value
SI	Sending current value (also instable)
SIR	Sending current value (also instable) and repeating
Т	Taring
Z	Zeroing

Extract of the mostly used KCP orders:

Example:

Order	S	
Possible replies	S_S100.00_ g S_I S_+ or S	Order accepted, execution of the order started, currently another order is executed, timeout reached, over- or underload

15.2 KERN alibi memory

For weighings where verification is mandatory and which are to be analysed and processed by a PC (e.g. printing out a packing list using a PC instead of a printer connected directly to the balance) electronic archiving is required by the metrological authorities by a verifiable data memory which cannot be manipulated. These stored data strings can be retrieved & displayed at any time via a connected PC

- The Alibi memory offers the possibility to store up to 250.000 weighing results, when the memory is exhausted, already used IDs are overwritten (starting with the first ID).
- By pressing the Print key or by KCP remote control command "S" or "MEMPRT" the storage process can be performed.
- The weight value (N, G, T), date and time and a unique alibi ID are stored.
- When using a print option, the unique alibi ID is also printed for identification purposes as well.
- The stored data can be retrieved via the KCP command "MEMQID". This can be used to query a specific single ID or a series of IDs.
- Example:
 - $\circ~$ MEMQID 15 \rightarrow The data record which is stored under ID 15 is returned.
 - $\circ~$ MEMQID 15 20 \rightarrow All data sets, which are stored from ID 15 to ID 20, are returned.

The detailed description can be found in the *KERN Communication Protocol* manual, available in the Downloads tab on the home page of KERN (<u>www.kern-sohn.com</u>).

1	 Protection of stored legally relevant data: After a record is stored, it will be read back immediately and be verified byte by byte. If error is found that record will be marked as an invalid record. If no error, then the record can be printed if needed. There is checksum protection stored in every record. All information on a printout is read from the memory with checksum verification, instead of direct from buffer.
	 Data loss prevention measures: The memory is write-disabled upon power-up. A write enable procedure is performed before writing a record to the memory. After a record is stored, a write disable procedure will be performed immediately (before verification).
	 The memory has a data retention period longer than 20 years
15.3 Issue functions

15.3.1 Add-up mode < └u미 >

With this function the individual weighing values are added into the summation memory by pressing a button and edited when an optional printer is connected.

Activate function:

1

- Solution Select the setting < □□> and confirm on → button.
- \Rightarrow To exit the menu, press the navigation key \leftarrow repeatedly

Condition: Menu setting

Add-up weighed goods:

- ⇒ If required, place empty container on scale and tare.
- ⇒ Place first good to be weighed on balance. Wait until stability display (▲ ▲) appears and then press the PRINT-button. The display changes to < □□□□ ↓>, followed by the current weighing value. The weighing value is stored and edited by the printer. The symbol ∑ pops up. Remove the weighed good.
- Place second good to be weighed on balance. Wait until stability display (►▲) appears and then press the PRINT-button. The display changes to < └└│2 >, followed by the current weighing value. The weighing value is stored and edited by the printer. Remove the weighed good.
- \Rightarrow Add-up more weighed goods as described above.
- \Rightarrow You can repeat this process until the capacity of the scales is exhausted.

Display and edit sum "Total":

⇒ Press the PRINT key long time. The number of weighings and the total weight are edited.

The sum memory is deleted; the symbol [. Σ .] extinguishes.

Sample log (KERN YKB-01N):

 $Menu \ setting < Pr \ \Pi adE > \rightarrow < Far \ \Pi AE > \rightarrow < Share >$

Sample log (KERN YKB-01N):

Menu setting

<PrNodE> > < UE (Ght> > < SGLPrt> > <on>

15.3.2 Data output after pressing the PRINT button < 미유고마유니 > Activate function:

- ⇒ In Setup menu invoke the menu setting < Pr in $E > \rightarrow < Pr \square dE > \rightarrow < Er$ i G > and confirm with → button.
- ⇒ For a manual data output select the menu setting < \\A□uAL > with the navigation keys ↓1 and confirm on the → button.
- \Rightarrow Use the navigation keys $\downarrow\uparrow$ to select the setting < $\Box \neg$ > and confirm on \rightarrow button.
- \Rightarrow To exit the menu, press the navigation key \leftarrow repeatedly.

Place goods to be weighed on balance:

- \Rightarrow If required, place empty container on scale and tare.
- ⇒ Place goods to be weighed. The weighing value is edited by pressing the PRINTbutton.

15.3.3 Automatic data output < Ruto >

Data output happens automatically without having to press the **PRINT** button as soon as the corresponding output condition has been met, dependent on the setting in the menu.

Enable function and set the output condition:

- ⇒ In Setup menu invoke the menu setting < Pr in $E > \rightarrow < Pr$ ∩ $\Box dE > \rightarrow < Er$ i $\Box >$ and confirm with \rightarrow button.
- For an automatic data output select the menu setting $< \exists u \exists u =$ using the navigation keys \downarrow 1 and confirm by the \rightarrow button.
- Use the navigation keys 11 to select the setting < $\Box \Box >$ and confirm on \rightarrow button. < $\Box \Box \Box \Box \Box \Box =$ is displayed.
- \Rightarrow Acknowledge by \rightarrow button and set the required output condition with the navigation keys \downarrow 1.
- \Rightarrow Acknowledge by \rightarrow button.
- \Rightarrow To exit the menu press the navigation key \leftarrow repeatedly.

Place goods to be weighed on balance:

- ⇒ If required, place empty container on scale and tare.
- ⇒ Place weighed goods and wait until the stability display (► →) appears. The weighing value is issued automatically.

15.3.4 Continuous data output < こロロヒ >

Enable function and set the output interval:

- ⇒ In Setup menu invoke the menu setting < Pr $nE > \rightarrow < Pr nDdE > \rightarrow < Er$ nE > and confirm with → button.
- For a continuous data output select the menu setting $< \Box \Box \Box \Box L >$ using the navigation keys \downarrow 1 and confirm on \rightarrow button.
- Solution Select the setting < □ □ > and confirm on → button.
- \Rightarrow < \square PEEd> is displayed.
- Acknowledge with the → button and set the required time interval with the navigation keys ↓1 (numerical input see chap. 3.2.2)
- \Rightarrow Set the required output condition <2E a > & <5EBbLE >.
- \Rightarrow To exit the menu press the navigation key \leftarrow repeatedly.

Place goods to be weighed on balance

- \Rightarrow If required, place empty container on scale and tare.
- \Rightarrow Place goods to be weighed.
- \Rightarrow The weighing values are issued according to the defined interval.

Sample log (KERN YKB-01N):

_	S	D	1.9997	40
	ŝ	D	1.9999	ka
	ŝ	D	1.99999	ko
	ŝ	D	1.99999	ka
	ŝ	S	2.0000	ka
	ŝ	5	2.0000	ka
	ŝ	<u> </u>	2.0000	ka
	ŝ	ŝ	2.0000	ka
	ŝ	5	1 9999	ko
	ŝ	D	4 9999	ka
	e e	5	2 0002	kg
	2	D D	2 1400	
	2	D D	2.4109	kg ka
	ĩ	D	2 0006	kg ka
	-	5	2. 2220	
	2	5	2. 3336	Kg
	2	5	2. 9997	K g
	2	0	2.9997	kg
	5	5	2.9996	Kg
			2.9996	кg

15.4 Data format

- ⇒ In the setup menu call up the menu setting < Pr $nE > \rightarrow < Pr n dE > \rightarrow < Br n dE$
- ⇒ Use the navigation keys \downarrow 1 to select the menu setting < For \square HE > and confirm on → button.
- Solution ⇒ Use the navigation buttons ↓↑ to select the desired setting. Options:

< Standard measuring protocol

<Lonu > Detailed measuring protocol

- \Rightarrow Confirm setting with \rightarrow button.
- \Rightarrow To exit the menu press the navigation key \leftarrow repeatedly.

Sample log (KERN YKB-01N):

Forf	18E → Shor	۰E	ForNAt - LonG
N: T: G:	SS	2.0000 kg 0.5000 kg 2.5000 kg	N: 5 D 2.0000 kg Tara weight after x: 0.5000 kg Gross weight: 2.5000 kg

16 Servicing, maintenance, disposal

Before any maintenance, cleaning and repair work disconnect the appliance from the operating voltage.

16.1 Cleaning

Please do not use aggressive cleaning agents (solvents or similar agents), but a cloth dampened with mild soap suds. Ensure that no liquid penetrates into the device. Polish with a dry soft cloth.

Loose residue sample/powder can be removed carefully with a brush or manual vacuum cleaner.

Spilled weighing goods must be removed immediately.

16.2 Servicing, maintenance

- ⇒ The appliance may only be opened by trained service technicians who are authorized by KERN.
- ⇒ Before opening, disconnect from power supply.

16.3 Disposal

Disposal of packaging and appliance must be carried out by operator according to valid national or regional law of the location where the appliance is used.

17 Instant help for troubleshooting

In case of an error in the program process, briefly turn off the balance and disconnect from power supply. The weighing process must then be restarted from the beginning.

Fault	Possible cause
The weight display does not glow.	The balance is not switched on.The mains supply connection has been interrupted
	(mains cable not plugged in/faulty).
	Power supply interrupted.
The displayed weight is per- manently changing	Draught/air movement
	Table/floor vibrations
	Weighing plate has contact with foreign objects.
	 Electromagnetic fields / static charging (choose dif- ferent location/switch off interfering device if possi- ble)
The weighing result is obvi-	 The display of the balance is not at zero
	 Adjustment is no longer correct.
	• The balance is on an uneven surface.
	Great fluctuations in temperature.
	Warm-up time was ignored.
	 Electromagnetic fields / static charging (choose dif- ferent location/switch off interfering device if possi-

ble)

18 Error messages

Error message	Explication	
ጋር ጊ ጉ	Zero setting range exceeded	
undErl	Zero setting range not achieved	
instAp	Load instable	
UronG	Adjustment error	
LJ	Underload	
[]	Overload	
LobAt	Capacity of batteries / rechargeable batteries exhausted	

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