

Measure Brix and Salt at the same time!

Save time and space.

It’s even equipped with a data communication function!

In just-a-minute

From 1 year to 2 years

Free Extended Warranty

1. It requires **only 1 minute!**

Simply by answering questions, warranty period is extended from 1 year to 2 years.


2. ATAGO Logger NFC can also be downloaded at the same time.



Access now ⇒

(The registration page can be accessed from ATAGO website.)


Parts



Preparation

Remove the tape from the battery compartment and close the cover.

Note The cover must be closed properly to maintain water resistance. Keep the batteries in during storage.



Measurement Principles

This instrument uses the electric conductivity method to measure and display salt concentrations % (g/100g). When complex samples containing ingredients other than salt are measured, the conductivity readings may be different from readings by other methods. Always dilute a complex sample to 10% by weight when its Brix exceeds 6%.

Brix is a measurement of the total dissolved solids (TDS) in a solution and measured by a refractometer. Check the Brix of you sample with a refractometer. For optimum results, it is recommended to dilute complex samples that are 6% Brix or higher.

Brix and Salt Measurements

Making Dilutions (Example of 10-fold dilution)


1.Place of sample.
Dilute by weight

10g of sample.

2. Add water

90g of water.

3. Gently stir.



Brix and Salt measurements

1. Clean the sensor.


2. Apply the sample.

At least 0.6mL

3. Measurements

Press the START button


Results



Offset feature use #1

Input a coefficient (a) of 10, and the value multiplied by 10 will be displayed.

Displayed value



Sample Preparation

Drinkable as is (less than 6% Brix)
⇒No dilution is necessary

Liquid condiments (over 6% Brix, over 10% salt, and high in non-salt components)
Soy sauce, Worcester sauce, etc.
⇒Please dilute.
See “Making Dilutions”

Paste
Mayonnaise, miso paste, ketchup, etc.
⇒Please dilute.
See “Making Dilutions”

Solid food
Pickles, ham, cheese, chips, etc.
⇒Please mince/grind and dilute.
See “Making Dilutions”

※Wait for approx. 5 minutes for the solids to settle to the bottom and measure the clear liquid on top.

Measurement Examples

Tomato puree 1.7%
Ketchup 3.0%
BBQ sauce 4.8%
Oyster sauce 9.4%
Salmon 2.4%
Salted cod roe 5.2%
Potage 1.2%
Miso soup 0.9%
Soy sauce 13.0%
Mayonnaise 1.6%
Pickled radish 3.6%
Pickles 1.7%
Ham 1.1%
Sausage 0.8%
Noodle soup 1.4%
Curry 1.6%
Gouda 0.9%
Butter 0.1%
Crackers 2.3%
Chips 1.4%

(Test data by ATAGO)

Preparation①Zero-setting

※Recommended on a daily basis.

Salt: Zero-setting

1. Clean the sensor (with nothing on the sensor).

2. Press the START button

3. Press the ZERO button

4. Apply tap water. At least 0.6mL

5. Press the START button.

6. Press the ZERO button.

Brix: Zero-setting

1. Clean the sensor (with nothing on the sensor).

2. Press the START button

3. Press the ZERO button

4. Apply tap water. At least 0.6mL

5. Press the START button.

6. Press the ZERO button.

Measurement : not 0.00%

Measurement : 0.00%

The instrument is zero-set correctly.

Proceed to measurement.

Salt: Offset Function

Discrepancies with titration

Due to the difference in measurement principles, readings from the instrument may not match up exactly with the readings by titration for certain samples. However, correlation between the two testing methods can be seen.


Offset feature use

Create a conversion chart between the two testing methods.

$y = a \times x + b$

y: titration readings
x: The instrument readings
a: coefficient (multiplication)
b: addition/subtraction number

An example of conversion chart



For addition/subtraction (b) [Range: -10.00 to 10.00]

1. Hold down ZERO while it is turned on. (At measured value displaying)

2. Press ZERO to select either addition (b) or subtraction (-b).

3. Press START to confirm.

4. Enter the addition/subtraction number.

5. Press START to confirm the addition/subtraction number.

“b” will appear

No plus sign will be displayed.

For addition, only “b” will appear.

The 1st digit: ZERO to change the number. START to confirm and move to the next digit.
The 2nd digit: ZERO to change the number. START to confirm and move to the next digit.
The 3rd digit: ZERO to change the number.

Next is to program a coefficient.

Coefficient

For coefficient (a) [Range: 0.50 to 5.00]

6. Enter the coefficient.

7. Press START to confirm the coefficient.

The 1st digit: ZERO to change the number. START to confirm and move to the next digit.
The 2nd digit: ZERO to change the number. START to confirm and move to the next digit.
The 3rd digit: ZERO to change the number.

The measurement range is shifted according to the offset settings.
Screen images when offset is on

During measurement (Salt)

After zero-setting with air

Offset “b” = addition of 0.30

For default Press the Heart button while setting up the Off-Set feature. Factory default value

Brix and Salt Measurements

Brix and Salt measurements

1. Clean the sensor.

2. Apply the sample.

At least 0.6mL

3. Measurements

Press the START button

Results

4. Low Sodium index

Press the Heart button

Low Sodium index

Cleaning

1. Wipe off the sample.

2. Clean with a mild soap, and rinse well with water.

3. Dry the area with tissues thoroughly.

Wash under running water (<50°C).

Do not submerge in water.

Do not use alcohol.

Error Messages

LLL

RRR

HHH

nnn

ooo

Err

Lo

- Brix : Sample not enough.
- Temperature : Below the range.
- Brix: Zero-setting with other than water.
- Salt: Calibration with other than reference solution.
- Above the measurement range.
- Brix : Too Bright. Shade the sample stage with your hand.
- Low sodium index error.
- Creating user scales error.
- The battery is low.

Salt: Checking with Reference Solution

When there is any doubt regarding accuracy of measurement results, adjust the reference value according to the following procedure.

memo The reference solution is available from ATAGO.

Part No. RE-120250 NaCl Solution 2.50%

Part No. RE-120900 NaCl Solution 9.00%

Part No. RE-121500 NaCl Solution 15.00%

memo Calibration automatically recognizes 2.50%, 9.00%, 15.00%.

[1]Salt : Zero-setting

1. Clean the sensor

2. Salt: Zero-setting

Press the ZERO button (with nothing on the sensor) to zero-set.

[2]Salt : Checking with Reference Solution (NaCl Solution2.50%, 9.00%, 15.00%)

1. Clean the sensor

2. Apply the reference solution.

At least 0.6mL

3. Press the START button

Press the START button

4. Results

Salt: 2.50±0.13
9.00±0.45
15.00±1.50

The instrument is calibrated correctly.

Within standard value range

Calibration is required.

[3]Salt : Calibration

1.Salt: Calibration

Calibration is complete

Apply the reference solution. Press the START and ZERO buttons simultaneously.

Inserting Batteries

1.Insert a coin in the groove on the battery compartment cover.
Turn the coin counterclockwise to remove the cover.

2.Insert batteries, observing the correct polarity.

3.Align the cover and push it down.

4.Close the battery compartment cover by pushing the cover in with a coin in the groove and turning it clockwise until it stops.
※Turning excessively may cause malfunction.

O-ring

NoteWhen the O-ring on the battery compartment cover is dirty or damaged, the water resistance may be compromised.

Safety Precautions

Read and follow all safety instructions before operating the instrument.

WARNING

- When measuring hazardous materials, use proper safety procedures, materials, and clothing to avoid personal injury. Anyone handling hazardous materials should understand its properties and its safety requirements.
- If the instrument is dropped or subjected to a strong impact, contact your supplier for inspection.
- Do not attempt to repair, modify, or disassemble the instrument.

CAUTION

- Before use, carefully read the instruction manual and fully understand the function and operation for each part of the instrument.
- ATAGO is not liable for any loss and damage caused by the measurement and use of this instrument.
- If this instrument is used to measure highly acidic samples, the sensor section and sample stage may be damaged, resulting in inaccurate measurements.
- Do not use any metal tools when applying sample to the sensor section. The metal can damage the sensor section. If the sensor section is scratched or damaged, inaccurate measurements will occur.
- When the unit needs to be washed, use water at a temperature not exceeding 50°C.
- Only use the specified battery type. Observe proper polarities, properly aligning the anodes and cathodes.
- Do not leave the instrument in a location exposed to direct sunlight or near a heat source for any extended period of time.
- Do not change the ambient temperature of the instrument suddenly.
- Do not place the instrument where it will be subject to strong vibrations.
- Do not use the instrument where there are excessive amounts of dust.
- Do not store the instrument in an extremely cool area.
- Do not set or drop heavy objects on top of the instrument.
- Loosen the battery compartment cover for air transportation.
- The instrument is water-resistant, not waterproof, and should not be submerged.

Contents

Main unit···1 Instruction Manual (this book)···1 AAA batteries···2

ATAGO instruments are rigorously inspected to ensure each unit meets the highest standards of quality assurance.

Automatic Temperature Compensation

The Automatic Temperature Compensation (ATC) feature is based on temperature detected by the thermo sensor located near the sensor area. ATC may not work correctly when the temperature of the sensor area is not the same as the actual temperature of the sample. When measuring a hot or cold sample, let it sit on the sensor for approximately 20 seconds and measure, or take multiple readings until measurements become stable.

Storage and Maintenance

Store the instrument in a dry place away from direct sunlight. Exposure to humidity may cause condensation inside, and exposure to direct sunlight may cause the plastic to warp.

- Cleaning** Clean and dry the sensor area thoroughly after use, leaving no sample residues or water.
(For oily samples:)
Remove oily residues with mild soap, and then, rinse with water.
- Storage** Store the instrument away from direct sunlight at a stable temperature with as little fluctuation as possible.

Repair and Warranty

The instrument is warranted for one year from the date of purchase. This warranty is void if the instrument shows evidence of the following. Send the included batteries as well if they are still in use.

- Having been disassembled by unauthorized personnel
- Damages to the sensor section
- Water damage or having been dropped
- Having been misused and/or operated outside the environmental specifications
- Leakage from batteries other than those included with the unit

Repair services are available for a fee after the warranty expires. Contact an ATAGO authorized service center for service and support.

Please have the serial number information ready when contacting a service center.

Specifications

Measurement range	Brix 0.0 to 90.0% , Temperature10.0 to 100°C SALT 0.00 to 15.00% (g/100g)
Resolution	Brix 0.1% SALT 0.01% Temperature 0.1°C
Measurement accuracy	Brix ±0.2% Temperature ±1°C SALT Displayed value ±0.05% (for salt concentration of 0.00 to 0.99%) Relative precision ±5% (for salt concentration of 1.00 to 9.99%) Relative precision ±10% (for salt concentration of 10.00 to 15.00%)
Measurement time	Approx. 3 seconds
Backlight	The backlight stays on for 30 seconds after any button is pressed.
Output	NFC Forum Type 4 Tag ISO/IEC 14443 Type A
Output category	Date Time, Brix[%],Salt[g/100g], Low Sodium Index, Temp[degC] (e.g.) 2021/01/17 09:30:45, 76.1,3.71,0.05,21.3
Maximum number of data history	100
Automatic temperature compensation range	Brix : 10 to 100°C Acid : 10 to 40°C
Ambient temperature range	10 to 40°C
International Protection class	IP65
Power supply	Two (2) AAA alkaline batteries
Dimensions and Weight	55(W)×31(D)×109(H)mm , 100g (main unit only)

About Data Transmission Function

This instrument stores maximum number of 100 measurement data. This instrument is equipped with NFC (Near Field Communication) technology. Data history can be accessed by bringing PAL-NFC to any Android devices, iPhone or PC-linked USB NFC Reader/Writer* (in conformance to PC/SC specification).
* Operation tested with SONY USB NFC Reader PaSoRi RC-S380.

Caution Data history exceeding 100 will overwrite old activity with new data, replacing the oldest recorded information first.

2017/08/17 09:30:45,12.3,20.4
2017/08/17 09:30:50,12.3,20.4
2017/08/17 09:30:55,12.4,20.4
...

Maximum of 100 measurement results are recorded with time stamps.

Example of data history

Android devices / iPhone

Laptop or PC + USB NFC Reader/Writer

Example of data history read out

Preparation

(1) Software installation Install a software to readout the NFC tag ahead of time.

Android devices / iPhone

Laptop or PC + USB NFC Reader/Writer

Applicable Application Software (app) "NFC Reader"

* If an NFC tag reader app is already installed on the Android devices or iPhone, this app can be used.

Data history can be exported to Microsoft^(R) Excel^(R)(for Windows^(R)) using a PAL NFC software "ATAGO Logger (NFC)."

* The app "ATAGO Logger (NFC)" is available for download : <http://www.atago.net/ur/>

2021/01/17 09:30:45	LLL	LLL	LLL	20.4
2021/01/17 09:31:50	12.5	0.31	0.39	20.5
2021/01/17 09:32:12	3.2	1.25	0.39	21.0
2021/01/17 09:34:26	AAA	AAA		21.2
2021/01/17 09:43:07	OOO			22.3
2021/01/17 09:43:18		OOO		22.5
2021/01/17 09:45:39	38.5	5.89	0.15	25.1
2021/01/17 09:46:07	39.2	5.92	0.15	25.3

LLL : Lower limit error.
HHH : Upper limit error.
OOO : Zero setting complete.
nnn : External light error.
AAA : Zero setting error.
ooo : Low Sodium index error.

A0135897
667937E581

Back side of the body
Instrument's serial number
Last 10 digits of NFC chip number (NFC'S serial number)
You can use the NFC'S serial number to identify which instrument the readings correspond to.

You can check the NFC chip number (serial number) by using an app that can read the serial number.
Note: Not all NFC apps capture the NFC'S serial number
Suggested app: "NFC Tools"

(2) Date and time setting

Set the date and time (year [the last two digits of the western calendar], month, date, time and minute) prior to data history readout.

memo Reset the date and time when batteries are removed for 24hours or more.

While powered on...
(Example image)

10's place

1's place

10's place

1's place

10's place

1's place

10's place

1's place

START

START

START

START

START

START

START

START

Confirm

Confirm

Confirm

Confirm

Confirm

Confirm

Confirm

Confirm

Year

Year

Month

Day

Day

Hour

Minute

Minute

0-9

0-9

01-12

0-3

0-9

00-23

0-5

0-9

End

Press the ZERO button until the desired number is displayed.

Pressing the ZERO button will change the number.

Top right screen display during date and time set up
(Seconds : Fixed 00)

Year : 99

Month : 12

Day : 31

Hour : 24

Minute : 00

Data history readout

Android devices / iPhone

Laptop or PC + USB NFC Reader/Writer

1. Launch NFC Reader (or other NFC tag reading app).
2. Position NFC on Android device or iPhone to the "NFC" logo at the bottom portion of the PAL then bring in contact.

* Do not move it.
(Hold for 1 second or more.)

1. Launch ATAGO Logger.
2. Bring the bottom part (where the "NFC" logo is) of PAL in contact with the NFC mark on the IC card reader/writer.

* Do not move it.
(Hold for 1 second or more.)

* Be sure to establish the PC and IC card reader/writer connection in advance by setting up (and installing the driver) IC card reader/writer.

* Data history can be read out by holding up the USB NFC Reader/writer to the PAL unit.

* NFC position on Android device or iPhone differs to the mode I.

All recorded data stored in this instrument are read out.

* If data history is not read out, bring both in contact and move the one that is over the other device in a forward and back or left and right in a small motion.

Delete data history

1. Quickly (3 seconds or less) do the following button operation.
(a) While pressing the START button, press the ZERO button two times.
(b) Quickly release the START button.

While powered on...
(Example image)

START button

START

ZERO button

ZERO

ZERO

0 second

1 second

2 seconds

3 seconds

START

START

START

Press the START button.

Press the ZERO button.

Delete data history

End without deleting the data history.

All data history will be deleted from this instrument.

Caution Deleted data history can not be restored.
memo A data history can not be selected.

ATAGO's instruments are rigorously inspected to ensure each unit meets the highest standards of quality assurance.

ATAGO CO.,LTD.

Headquarters: The Front Tower Shiba Koen, 23rd Floor
2-6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan
TEL: 81-3-3431-1943 overseas@atago.net <http://www.atago.net/>

ATAGO U.S.A., Inc.

TEL: 1-425-637-2107 customerservice@atago-usa.com

ATAGO INDIA Instruments Pvt.Ltd.

TEL: 91-22-28544915 / 40713232
customerservice@atago-india.com

ATAGO THAILAND Co.,Ltd.

TEL: 66-21948727-9 ,66-21171549
customerservice@atago-thailand.com

ATAGO BRASIL Ltda.

TEL: 55 16 3913-8400 customerservice@atago-brasil.com

ATAGO ITALIA s. r. l.

TEL: 39 02 36557267 customerservice@atago-italia.com

ATAGO CHINA Guangzhou Co.,Ltd.

TEL: 86-20-38108256 info@atago-china.com

ATAGO RUSSIA Ltd.

TEL: 7-812-777-96-96 info@atago-russia.com

ATAGO NIGERIA Scientific Co.,Ltd.

TEL: 234-707-558-1552 atagonigeria@atago.net

ATAGO KAZAKHSTAN Ltd.

TEL: 7-727-257-08-95 info@atago-kazakhstan.com
2102K Printed in Japan